

DNA Sample Preparation

Andrea Slade
Veeco Life Science Applications Scientist

1. DNA: Lambda *Bst*E II Digest (Sigma Aldrich Catalog # D-9793)
2. Storage Buffer: 10mM Tris-HCl (pH 8.0), 1 mM EDTA
3. Adsorption Buffer: 40 mM HEPES (pH 6.8), 10 mM NiCl₂

NOTE: Nickel Chloride (NiCl₂) can be a skin irritant in both solid and solution form. Please wear gloves when handling any NiCl₂ containing substances and wash hands when finished. Please, dispose of any NiCl₂ containing solutions in the appropriate labeled waste container – on lab bench in Biology Lab. Do NOT dispose in sink or in trash bins.

Storage Buffer **(10mM Tris-HCl (pH 8.0), 1 mM EDTA)**

10 mM Tris-HCl (MW = 157.60 g/mol)
1 mM EDTA (MW = 372.24 g/mol)
50 mL buffer volume

For Tris-HCl,
 $m = (10 \text{ mmol/L})(0.050 \text{ L})(157.60 \text{ mg/mL}) = 0.0788 \text{ g}$

For EDTA,
 $m = (1 \text{ mmol/L})(0.050 \text{ L})(372.24 \text{ g/mol}) = 0.0186 \text{ g}$

- a. Weigh out 0.0788 g of TRis-HCl and 0.0186 g of EDTA.
- b. Dissolve in 50 mL of Milli-Q water
- c. Adjust pH to 8.0 (using 1 M NaOH)
- d. Store buffer in fridge when not in use.
- e. Note: Buffer must be filtered (using 0.45µm syringe filter) before used to dilute DNA.

Adsorption Buffer **(40 mM HEPES (pH 6.8), 10 mM NiCl₂)**

40 mM HEPES (MW = 238.3 g/mol)
1 mM NiCl₂ (MW = 237.7 g/mol)
50 mL buffer volume

For HEPES,

$$m = (40 \text{ mmol/L})(0.050 \text{ L})(157.60 \text{ mg/mL}) = 0.4766 \text{ g}$$

For NiCl_2 ,

$$m = (10 \text{ mmol/L})(0.050 \text{ L})(372.24 \text{ g/mol}) = 0.1188 \text{ g}$$

- Weigh out 0.4766 g of HEPES and 0.1188 g of NiCl_2 .
- Dissolve in 50 mL of Milli-Q water
- Adjust pH to 6.8 (using 1 M NaOH)
- Store buffer in fridge when not in use.
- Note: Buffer must be filtered (using 0.45 μm syringe filter) before used to dilute DNA.

DNA Preparation

- DNA Solution comes as 863 $\mu\text{g/mL}$ solution (~400-500 μL volume)
- Mix solution by finger tapping the tube of solution (Do NOT vortex the solution.)
- Aliquot solution into eppendorf tubes as 10 μL / tube.
- Store the eppendorf tubes in the freezer.

DNA Dilution (863 $\mu\text{g/mL}$ \rightarrow 3 $\mu\text{g/mL}$):

- Take one 10 μL aliquot of stock DNA and thaw to room temperature.
- Remove Storage Buffer (10 mM Tris-HCl (pH 8.0), 1 mM EDTA) from fridge and allow to come to room temperature.
- Add 2.87 mL of Storage Buffer to the 10 μL of DNA for a final concentration of 3 $\mu\text{g/mL}$.

$$v = [(863 \mu\text{g/mL})(0.010 \text{ mL})] / (3 \mu\text{g/mL}) = 2.88 \text{ mL}$$
$$(10 \mu\text{L DNA}) + (2.87 \text{ mL Storage Buffer}) = 2.88 \text{ mL total volume}$$

- Mix the solution thoroughly by finger tapping the tube.
- Aliquot the 3 $\mu\text{g/mL}$ of DNA solution into eppendorf tubes as 100 μL / tube.
- Store the eppendorf tubes in the freezer until ready to use in experiment.

On Day of Experiment...

DNA Dilution (3 $\mu\text{g}/\text{mL}$ \rightarrow 0.5 $\mu\text{g}/\text{mL}$):

- a. Remove 1 eppendorf tube of 3 $\mu\text{g}/\text{mL}$ DNA solution from the freezer and thaw to room temperature.
- b. Add 500 μL of Adsorption Buffer (40 mM HEPES (pH 6.8), 10 mM NiCl_2) to the tube of 100 μL of DNA for a final concentration of 0.5 $\mu\text{g}/\text{mL}$. (Make sure buffer has been filtered: See buffer section, Step e)

$$v = [(3 \mu\text{g}/\text{mL})(0.100 \text{ mL})] / (0.5 \mu\text{g}/\text{mL}) = 600 \mu\text{L}$$
$$(100 \mu\text{L DNA}) + (500 \mu\text{L Adsorption Buffer}) = 600 \mu\text{L total volume}$$

- c. Mix the solution thoroughly by repeatedly pipetting the solution in and out of the same pipette tip.
- d. Pipette 100 μL of solution onto a freshly cleaved mica substrate (mica + Teflon disk glued to a magnetic puck/glass slide).
- e. Allow the DNA to adsorb onto the mica surface for 5 minutes.
- f. Rinse the mica with Adsorption Buffer to remove DNA that has not adsorbed to the surface.
- g. Add fresh Adsorption Buffer to form a sufficient size droplet of fluid and begin imaging.

This procedure should give you ~20-30% surface coverage of DNA over a 1 μm scan area.

Place remaining solution in fridge and use with 3 days. Discard after 3 days and use fresh DNA from freezer (3 $\mu\text{g}/\text{mL}$) and dilute.